





ZCCHC4 Antibody

Product Code	CSB-PA887984LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5U6
Immunogen	Recombinant Human Zinc finger CCHC domain-containing protein 4 protein (25-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger CCHC domain-containing protein 4, ZCCHC4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZCCHC4

Image



Western blot
 All lanes: ZCCHC4 antibody at 0.6µg/ml + A431 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 60, 26, 59 kDa
 Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA887984LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.