



ZDHHC5 Antibody, FITC conjugated

Product Code	CSB-PA861196LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0B5
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC5 protein (231-500AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC5 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 5) (DHHC-5) (Zinc finger protein 375), ZDHHC5, KIAA1748 ZNF375
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZDHHC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.