



ZG16 Antibody, Biotin conjugated

Product Code	CSB-PA528909LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60844
Immunogen	Recombinant Human Zymogen granule membrane protein 16 protein (17-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zymogen granule membrane protein 16 (Zymogen granule protein 16) (hZG16) (Secretory lectin ZG16), ZG16
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZG16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.