



ZNF232 Antibody, HRP conjugated

Product Code	CSB-PA026599LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNY5
Immunogen	Recombinant Human Zinc finger protein 232 protein (1-417AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 232 (Zinc finger and SCAN domain-containing protein 11), ZNF232, ZSCAN11
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF232
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.