





ZNF490 Antibody

Product Code	CSB-PA026786LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULM2
Immunogen	Recombinant Human Zinc finger protein 490 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 490, ZNF490, KIAA1198
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF490

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA026786LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA026786LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.