




Il18bp Antibody

Product Code	CSB-PA16769A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Z0M9
Immunogen	Recombinant Mouse Interleukin-18-binding protein (29-193AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-18-binding protein (IL-18BP) (Interferon gamma-inducing factor-binding protein), Il18bp, Igifbp
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Il18bp

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA16769A0Rb at dilution of 1:20

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.