



SLC39A1 Antibody, FITC conjugated

Product Code	CSB-PA17954C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY26
Immunogen	Recombinant Human Zinc transporter ZIP1 protein (126-179AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter ZIP1 (Solute carrier family 39 member 1) (Zinc-iron-regulated transporter-like) (Zrt- and Irt-like protein 1) (ZIP-1) (hZIP1), SLC39A1, IRT1 ZIP1 ZIRTL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC39A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.