



TNFRSF10B Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA08165C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O14763 |
| Immunogen | Recombinant Human Tumor necrosis factor receptor superfamily member 10B protein (56-210AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Tumor necrosis factor receptor superfamily member 10B (Death receptor 5) (TNF-related apoptosis-inducing ligand receptor 2) (TRAIL receptor 2) (TRAIL-R2) (CD antigen CD262), TNFRSF10B, DR5 KILLER TRAILR2 TRICK2 ZTNFR9 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | TNFRSF10B |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |