



TNFRSF10B Antibody, Biotin conjugated

Product Code	CSB-PA08165D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14763
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 10B protein (56-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 10B (Death receptor 5) (TNF-related apoptosis-inducing ligand receptor 2) (TRAIL receptor 2) (TRAIL-R2) (CD antigen CD262), TNFRSF10B, DR5 KILLER TRAILR2 TRICK2 ZTNFR9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFRSF10B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.