





Recombinant *Macaca fascicularis* Transthyretin (TTR)

Product Code	CSB-MP025270MOV
Relevance	Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain (By similarity).
Abbreviation	Recombinant Cynomolgus monkey TTR protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8HXL1
Product Type	Recombinant Proteins
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (Cynomolgus monkey)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GPTGVDESKCPLMKVLDVAVRGSPAVNVAVNVFKKAADETWAPFASGKTSES GELHGLTTEEEFVEGIYKVEIDTKSYWKSLSGSPFHEHAEVVFTANDSGPRHYT IAALLSPYSYSTTAVVTNPKE
Research Area	Cardiovascular
Source	Mammalian cell
Target Names	TTR
Expression Region	21-147aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Myc-tagged
Mol. Weight	17.7kDa
Protein Length	Full Length of Mature Protein


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of Mammalian Cell host and target protein, the LC-MS/MS Analysis result of CSB-MP025270MOV could indicate that this peptide derived from Mammalian Cell-expressed Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey) TTR.



Based on the SEQUEST from database of Mammalian Cell host and target protein, the LC-MS/MS Analysis result of CSB-MP025270MOV could indicate that this peptide derived from Mammalian Cell-expressed Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey) TTR.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.