




Recombinant Mouse Low affinity immunoglobulin gamma Fc region receptor III (Fcgr3), partial

Product Code	CSB-MP357412MO
Relevance	Receptor for the Fc region of complexed immunoglobulins gamma. Low affinity receptor.
Abbreviation	Recombinant Mouse Fcgr3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08508
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ALPKAVVKLDPPWIQVLKEDMVTLMCEGTHNPGNSSTQWFHNGRSIRSQVQA SYTFKATVNDSGEYRCQMEQTRLSDPVDLGVISDWLLQTPQRVFLEGETITL RCHSWRNKLLNRISFFHNEKSVRYHHYKSNFSIPKANHSHSGDYCKGSLGS TQHQSKPVTITVQDPATTSSISLVWYHT
Research Area	Others
Source	Mammalian cell
Target Names	Fcgr3
Expression Region	31-215aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Myc-tagged
Mol. Weight	25.2kDa
Protein Length	Extracellular Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.