




Recombinant Human Cyclin-A1 Isoform 2 (CCNA1), partial

Product Code	CSB-BP004804HU
Relevance	May be involved in the control of the cell cycle at the G1/S (start) and G2/M (mitosis) transitions. May primarily function in the control of the germline meiotic cell cycle and additionally in the control of mitotic cell cycle in some somatic cells.
Abbreviation	Recombinant Human CCNA1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P78396
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	METGFPAIMYPGSGFIGGWGEEYLSWEGPGLPDFVFQQQPVESEAMHCSNPK SGVVLATVARGPDACQILTRAPLGQDPPQRTVLGLLTANGQYRRTCQGGITRI RCYSGSENAFPPAGKKALPDCGVQEPKQGFDIYMDELEQGDSDSCSVREG MAFEDVYEVDTGTLKSDLHFLDFNTVSPMLVDSSLLSQSEDISSLGTDVINVT EYAEIYQYLREAEIRHRPKAHYMKKQPDITEGMRTILVDWLVEVGEEYKLRAE TLYLAVNFLDRFLSCMSVLRGKLQLVGTAAMLLASKYEEIYPPEVDEFVYITDD TYTKRQLLMEHLLLKVLAFDLTVPTTNQFLLQYLRRQGVCVRTENLAKYVAEL SLLEADPFLKYLPSLIAAAAFCLANYTVNKHFWPETLAAFTGYSLSSEIVPCLSEL HKAYLDIPHRPQQAIREKYKASKYLCVSLMEPPAVLLL
Research Area	Cancer
Source	Baculovirus
Target Names	CCNA1
Expression Region	1-464aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	54.2kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.