




Recombinant Human Signal-regulatory protein beta-1 (SIRPB1), partial

Product Code	CSB-MP021335HU
Relevance	Immunoglobulin-like cell surface receptor involved in the negative regulation of receptor tyrosine kinase-coupled signaling processes. Participates also in the recruitment of tyrosine kinase SYK.
Abbreviation	Recombinant Human SIRPB1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O00241
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	EDELQVIQPEKSVSVAAGESATLRCAMTSLIPVGPIMWFRGAGAGRELIYNQK EGHFPRVTTVSELTKRNNLDFSISISNITPADAGTYCYVKFRKGGSPDDVEFKSG AGTELSVRAKPSAPVVSGPAVRATPEHTVSFTCESHGFSPRDITLKWFKNGNE LSDFQTNVDPAGDSVSYSIHSTARVVLTRGDVHSQVICEIAHITLQGDLRGTA NLSEAIRVPPTLEVTQQPMRAENQANVTCQVSNFYPRGLQLTWLENGNVSRT ETASTLIENKDGTYNWMWLLVNTCAHRDDVLTQVEHDGQQAVSKSYALE ISAHQKEHGSDITHEAALAPTAPL
Research Area	Signal Transduction
Source	Mammalian cell
Target Names	SIRPB1
Expression Region	30-371aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Myc-tagged
Mol. Weight	41.3kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.