




# Recombinant Rat Tryptophan 2,3-dioxygenase UniRule annotation (Tdo2)

<b>Product Code</b>	CSB-EP023351RA
<b>Relevance</b>	Incorporates oxygen into the indole moiety of tryptophan. Has a broad specificity towards tryptamine and derivatives including D- and L-tryptophan, 5-hydroxytryptophan and serotonin.
<b>Abbreviation</b>	Recombinant Rat Tdo2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P21643
<b>Alias</b>	Tryptamin 2,3-dioxygenase Tryptophan oxygenase Tryptophan pyrrolase Tryptophanase
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MSGCPFSGNSVG YTLKNLSMEDNEEDGAQTGVNRASKGGLIYGDY LQLEKILNAQELQSEIKGNKIHDEHLFIITHQAYELWFKQILWELDSVREIFQNGHVRDERNMLKVMTRMHRVVVIFKLLVQQFSVLETMTALDFNDFREYLSPASGFQSLQFR LLENKIGVLQSLRVPYNRKH YRDNFEGDYNELLLKSEQEQ TLLQLVEAWLERTPGLEPHGFNFWGKFEKNILKGLEEEFLKIQAKKDSEEKEEQMAEFRKQKEVLLCLFDEKRHDYLLSKGERRLSYRALQGALMIYFYREEPRFQVPFQLLTSLMDIDT LMTKWRYNHVCMVHRMLGSKAGTGGSSGYYYLRSTVSDRYKVFVDFLNLSS YLVPRHWIPKMNPIIHKFLYTAEYSDSSYFSSDES D
<b>Research Area</b>	Metabolism
<b>Source</b>	E.coli
<b>Target Names</b>	Tdo2
<b>Expression Region</b>	1-406aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	63.9kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.