




# Recombinant Influenza A virus Hemagglutinin (HA), partial

<b>Product Code</b>	CSB-EP607773IER
<b>Relevance</b>	Binds to sialic acid-containing receptors on the cell surface, bringing about the attachment of the virus particle to the cell. This attachment induces virion internalization of about two third of the virus particles through clathrin-dependent endocytosis and about one third through a clathrin- and caveolin-independent pathway. Plays a major role in the determination of host range restriction and virulence. Class I viral fusion protein. Responsible for penetration of the virus into the cell cytoplasm by mediating the fusion of the membrane of the endocytosed virus particle with the endosomal membrane. Low pH in endosomes induces an irreversible conformational change in HA2, releasing the fusion hydrophobic peptide. Several trimers are required to form a competent fusion pore
<b>Abbreviation</b>	Recombinant Influenza A virus Hemagglutinin protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0HD60
<b>Alias</b>	Cleaved into the following 2 chains: Hemagglutinin HA1 chain Hemagglutinin HA2 chain
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Influenza A virus (strain A/Hickox/1940 H1N1)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	DTICIGYHANNSTDTVDTVLEKNVTVTHSVNLLLED SHNGKLCRLKGIAPLQLGK CNIAGWILGNPECESLLSKRSWSYIAETPNSENGTCYPGDFADYEELREQLSS VSSFERFEIFPKERSWPNHNINIGVTAACSHAGKSSFYKNLLWLTEKDGSYPNL NKS YV NK KE KE VLVLWGVHHP SNIENQKTLYRKENAYVSVVSSNYNRRFTPEI AERPKVRGQAGRMNYYWTLLEPGDTIIFEANGNLIAPWYAFALSRGLGSGIITS NASMDECDTKCQTPQGAINSSLPFQNIHPFTIGCEPKYVRSTKLRMVTGLRNIP SIQS
<b>Source</b>	E.coli
<b>Target Names</b>	HA
<b>Expression Region</b>	18-343aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	52.6kDa

**Protein Length**

Partial

**Image**

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.