




# Recombinant Triticum aestivum Non-specific lipid-transfer protein

<b>Product Code</b>	CSB-EP2021TQN
<b>Relevance</b>	Plant non-specific lipid-transfer proteins transfer phospholipids as well as galactolipids across membranes. May play a role in wax or cutin deposition in the cell walls of expanding epidermal cells and certain secretory tissues.
<b>Abbreviation</b>	Recombinant Triticum aestivum Non-specific lipid-transfer protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P24296
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Triticum aestivum
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	DCGHVDSLVRPCLSYVQGGPGPSGQCCDGVKNLHNQARSQSDRQSACNCL KGIARGIHNLNEDNARSIPPKCGVNLPLYTISLNIDCSR
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	N/A
<b>Expression Region</b>	25-113aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	25.5kDa
<b>Protein Length</b>	Full Length of Mature Protein

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.