




Recombinant *Aspergillus fumigatus* Peroxiredoxin Asp f3 (aspf3)

Product Code	CSB-EP525004NGS
Abbreviation	Recombinant <i>Neosartorya fumigata</i> aspf3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43099
Product Type	Recombinant Proteins
Immunogen Species	<i>Neosartorya fumigata</i> (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MSGLKAGDSFPSDVVFSYIPWSEDKGEITACGIPINYNASKEWADKKVILFALP GAFTPVCSARHVPEYIEKLPEIRAKGVDVVAVLAYNDAYVMSAWGKANQVTG DDILFLSDPDARFSK SIGWADEEGRTKRYALVIDHGKITYAALEPAKNHLEFSSA ETVLKHL
Research Area	Allergen
Source	<i>E.coli</i>
Target Names	aspf3
Expression Region	1-168aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	34.5kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.