



MDR1 Antibody

Product Code	CSB-PA001046HA01EKM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16875
Immunogen	Recombinant Entamoeba histolytica Multidrug resistance protein 1 protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Entamoeba histolytica
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multidrug resistance protein 1 (P-glycoprotein) (Fragment), MDR1
Immunogen Species	Entamoeba histolytica
Research Area	Others
Target Names	MDR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.