



ACOX3 Antibody

Product Code	CSB-PA001174LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15254
Immunogen	Recombinant Human Peroxisomal acyl-coenzyme A oxidase 3 protein (401-624AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal acyl-coenzyme A oxidase 3 (EC 1.3.3.6) (Branched-chain acyl-CoA oxidase) (BRCACox) (Pristanoyl-CoA oxidase), ACOX3, BRcox, PRcox
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ACOX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.