




ADAMDEC1 Antibody

Product Code	CSB-PA001297LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15204
Immunogen	Recombinant Human ADAM DEC1 protein (206-470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADAM DEC1 (EC 3.4.24.-) (A disintegrin and metalloproteinase domain-like protein decysin-1) (ADAM-like protein decysin-1), ADAMDEC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ADAMDEC1

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA001297LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.