





AP3M2 Antibody

Product Code	CSB-PA001878LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53677
Immunogen	Recombinant Human AP-3 complex subunit mu-2 protein (1-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AP-3 complex subunit mu-2 (Adaptor-related protein complex 3 subunit mu-2) (Clathrin assembly protein assembly protein complex 3 mu-2 medium chain) (Clathrin coat assembly protein AP47 homolog 2) (Clathrin coat-associated protein AP47 homolog 2) (Golgi adaptor AP-1 47 kDa protein homolog 2) (HA1 47 kDa subunit homolog 2) (Mu3B-adaptin) (P47B), AP3M2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AP3M2

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA001878LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA001878LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.