



APOOL Antibody, Biotin conjugated

Product Code	CSB-PA001949LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXV4
Immunogen	Recombinant Human MICOS complex subunit MIC27 protein (163-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MICOS complex subunit MIC27 (Apolipoprotein O-like) (Protein FAM121A), APOOL, CXorf33 FAM121A MIC27
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	APOOL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.