



ARRB1 Antibody, HRP conjugated

Product Code	CSB-PA002134LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49407
Immunogen	Recombinant Human Beta-arrestin-1 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-arrestin-1 (Arrestin beta-1) (Non-visual arrestin-2), ARRB1, ARR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARRB1