



ASGR2 Antibody, Biotin conjugated

Product Code	CSB-PA002208LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07307
Immunogen	Recombinant Human Asialoglycoprotein receptor 2 protein (63-286AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Asialoglycoprotein receptor 2 (ASGP-R 2) (ASGPR 2) (C-type lectin domain family 4 member H2) (Hepatic lectin H2) (HL-2), ASGR2, CLEC4H2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ASGR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.