





BCAT1 Antibody

Product Code	CSB-PA002600EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54687
Immunogen	Recombinant Human Branched-chain-amino-acid aminotransferase, cytosolic protein (4-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Branched-chain-amino-acid aminotransferase, cytosolic (BCAT(c)) (EC 2.6.1.42) (Protein ECA39), BCAT1, BCT1 ECA39
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BCAT1


Image



Immunohistochemistry analysis of human melanoma using CSB-PA002600EA01HU at dilution of 1:100



Immunohistochemistry analysis of human bladder cancer using CSB-PA002600EA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA002600EA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.