





BGN Antibody

Product Code	CSB-PA002683EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21810
Immunogen	Recombinant Human Biglycan protein (38-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Biglycan (Bone/cartilage proteoglycan I) (PG-S1), BGN, SLRR1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BGN


Image



Immunohistochemistry analysis of human pancreatic cancer using CSB-PA002683EA01HU at dilution of 1:100



Immunohistochemistry analysis of human placenta tissue using CSB-PA002683EA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA002683EA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)