




LINC01553 Antibody

Product Code	CSB-PA002950LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4QN01
Immunogen	Recombinant Human Putative uncharacterized protein encoded by LINC01553 protein (1-128AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative uncharacterized protein encoded by LINC01553, LINC01553, C10orf40
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LINC01553

Image



Western blot
 All lanes: LINC01553 antibody at 12µg/ml
 Lane 1: Rat heart tissue
 Lane 2: Mouse skeletal muscle tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 16 kDa
 Observed band size: 16 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.