




CFAP410 Antibody

Product Code	CSB-PA003758LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43822
Immunogen	Recombinant Human Cilia- and flagella-associated protein 410 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cilia- and flagella-associated protein 410 (C21orf-HUMF09G8.5) (Leucine-rich repeat-containing protein 76) (YF5/A2), CFAP410, C21orf2, LRRC76
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CFAP410

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA003758LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.