



C4BPA Antibody, HRP conjugated

Product Code	CSB-PA003951LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04003
Immunogen	Recombinant Human C4b-binding protein alpha chain protein (49-597AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C4b-binding protein alpha chain (C4bp) (Proline-rich protein) (PRP), C4BPA, C4BP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C4BPA