



CD300LB Antibody

Product Code	CSB-PA004919LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8K4G0
Immunogen	Recombinant Human CMRF35-like molecule 7 protein (18-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CMRF35-like molecule 7 (CLM-7) (CD300 antigen-like family member B) (CMRF35-A2) (Immune receptor expressed on myeloid cells 3) (IREM-3) (Leukocyte mono-Ig-like receptor 5) (Triggering receptor expressed on myeloid cells 5) (TREM-5) (CD antigen CD300b), CD300LB, CD300B CLM7 CMRF35A2 IREM3 LMIR5 TREM5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CD300LB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.