




CD8A Antibody

Product Code	CSB-PA004966LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01732
Immunogen	Recombinant Human T-cell surface glycoprotein CD8 alpha chain protein (22-182AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD8 alpha chain (T-lymphocyte differentiation antigen T8/Leu-2) (CD antigen CD8a), CD8A, MAL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD8A

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA004966LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.