




CDK2AP2 Antibody

Product Code	CSB-PA005063LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75956
Immunogen	Recombinant Human Cyclin-dependent kinase 2-associated protein 2 protein (1-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 2-associated protein 2 (CDK2-associated protein 2) (DOC-1-related protein) (DOC-1R), CDK2AP2, DOC1R
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDK2AP2

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA005063LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.