




CEACAM7 Antibody

Product Code	CSB-PA005167LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14002
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 7 protein (36-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carcinoembryonic antigen-related cell adhesion molecule 7 (Carcinoembryonic antigen CGM2), CEACAM7, CGM2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CEACAM7

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA005167LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.