



CORO7 Antibody

Product Code	CSB-PA005819LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57737
Immunogen	Recombinant Human Coronin-7 protein (630-900AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coronin-7 (Crn7) (70 kDa WD repeat tumor rejection antigen homolog), CORO7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CORO7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.