



CORO7 Antibody

Product Code	CSB-PA005819LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57737
Immunogen	Recombinant Human Coronin-7 protein (630-900AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coronin-7 (Crn7) (70 kDa WD repeat tumor rejection antigen homolog), CORO7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CORO7