



Cr2 Antibody

Product Code	CSB-PA005934LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19070
Immunogen	Recombinant Mouse Complement receptor type 2 protein (729-963AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement receptor type 2 (Cr2) (Complement C3d receptor) (CD antigen CD21), Cr2
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cr2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.