





CRABP2 Antibody

Product Code	CSB-PA005936LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29373
Immunogen	Recombinant Human Cellular retinoic acid-binding protein 2 protein (1-138AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cellular retinoic acid-binding protein 2 (Cellular retinoic acid-binding protein II) (CRABP-II), CRABP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CRABP2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA005936LA01HU at dilution of 1:100

**Western Blot**

Positive WB detected in: MCF-7 whole cell lysate

All lanes: CRABP2 antibody at 3 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 16 kDa

Observed band size: 16 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.