



CYSLTR1 Antibody, HRP conjugated

Product Code	CSB-PA006465LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y271
Immunogen	Recombinant Human Cysteinyl leukotriene receptor 1 protein (298-337AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cysteinyl leukotriene receptor 1 (CysLTR1) (Cysteinyl leukotriene D4 receptor) (LTD4 receptor) (G-protein coupled receptor HG55) (HMTMF81), CYSLTR1, CYSLT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CYSLTR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.