




BRINP1 Antibody

Product Code	CSB-PA006515LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60477
Immunogen	Recombinant Human BMP/retinoic acid-inducible neural-specific protein 1 protein (20-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BMP/retinoic acid-inducible neural-specific protein 1 (Deleted in bladder cancer protein 1), BRINP1, DBC1 DBCCR1 FAM5A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BRINP1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA006515LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA006515LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)