



DCN Antibody, FITC conjugated

Product Code	CSB-PA006554LC01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21793
Immunogen	Recombinant Bovine Decorin protein (31-360AA)
Raised In	Rabbit
Species Reactivity	Bovine
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Decorin (Bone proteoglycan II) (PG-S2), DCN
Immunogen Species	Bos taurus (Bovine)
Research Area	Signal Transduction
Target Names	DCN