



Defa5 Antibody

Product Code	CSB-PA006657EA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28312
Immunogen	Recombinant Mouse Alpha-defensin 5 protein (59-93AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-defensin 5 (Defensin-related cryptdin-5), Defa5, Defcr5
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Defa5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.