




DGKG Antibody

Product Code	CSB-PA006836LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49619
Immunogen	Recombinant Human Diacylglycerol kinase gamma protein (1-255AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diacylglycerol kinase gamma (DAG kinase gamma) (EC 2.7.1.107) (Diglyceride kinase gamma) (DGK-gamma), DGKG, DAGK3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DGKG

Image



Western blot
 All lanes: DGKG antibody at 4µg/ml + Mouse brain tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 90, 87, 85 kDa
 Observed band size: 90 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.