





DNAJC6 Antibody

Product Code	CSB-PA007040LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75061
Immunogen	Recombinant Human Putative tyrosine-protein phosphatase auxilin protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative tyrosine-protein phosphatase auxilin (EC 3.1.3.48) (DnaJ homolog subfamily C member 6), DNAJC6, KIAA0473
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DNAJC6


Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA007040LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA007040LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA007040LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)