



DNALI1 Antibody, FITC conjugated

Product Code	CSB-PA007047LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14645
Immunogen	Recombinant Human Axonemal dynein light intermediate polypeptide 1 protein (1-258AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Axonemal dynein light intermediate polypeptide 1 (Inner dynein arm light chain, axonemal) (hp28), DNALI1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNALI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.