



DNM1L Antibody, HRP conjugated

Product Code	CSB-PA007078EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00429
Immunogen	Recombinant Human Dynamin-1-like protein (496-710AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynamin-1-like protein (EC 3.6.5.5) (Dnm1p/Vps1p-like protein) (DVLp) (Dynamin family member proline-rich carboxyl-terminal domain less) (Dymple) (Dynamin-like protein) (Dynamin-like protein 4) (Dynamin-like protein IV) (HdynIV) (Dynamin-related protein 1), DNM1L, DLP1, DRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DNM1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.