



DSCR4 Antibody

Product Code	CSB-PA007197LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56555
Immunogen	Recombinant Human Down syndrome critical region protein 4 protein (1-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Down syndrome critical region protein 4 (Down syndrome critical region protein B), DSCR4, DCRB DSCRB
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DSCR4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.