



# DUS4L Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA007237LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O95620
<b>Immunogen</b>	Recombinant Human tRNA-dihydrouridine(20a/20b) synthase [NAD(P)+]-like protein (1-196AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	tRNA-dihydrouridine(20a/20b) synthase [NAD(P)+]-like (EC 1.3.1.-) (pp35) (tRNA-dihydrouridine synthase 4-like), DUS4L
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	DUS4L
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.