





EIF6 Antibody

Product Code	CSB-PA007582LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56537
Immunogen	Recombinant Human Eukaryotic translation initiation factor 6 protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 6 (eIF-6) (B(2)GCN homolog) (B4 integrin interactor) (CAB) (p27(BBP)), EIF6, EIF3A ITGB4BP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EIF6


Image




Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA007582LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA007582LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA007582LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry analysis of human small intestine tissue using CSB-PA007582LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.