



EMP1 Antibody, FITC conjugated

Product Code	CSB-PA007648LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54849
Immunogen	Recombinant Human Epithelial membrane protein 1 protein (22-66AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epithelial membrane protein 1 (EMP-1) (CL-20) (Protein B4B) (Tumor-associated membrane protein), EMP1, B4B TMP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EMP1