



EMP3 Antibody, FITC conjugated

Product Code	CSB-PA007650LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54852
Immunogen	Recombinant Human Epithelial membrane protein 3 protein (25-65AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epithelial membrane protein 3 (EMP-3) (Hematopoietic neural membrane protein 1) (HNMP-1) (Protein YMP), EMP3, YMP
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	EMP3